



## FEATURES

- Fabric Material: 65% Polyester / 35% Cotton
- Reusable, washable material
- Open back with ties provides easy don and doff
- Kimono style gown features curved armholes with elastic on sleeve ends, containing natural rubber latex
- Size: OSFM (one size fits most)
- Fabric made in the USA†

Reusable, isolation gown acts as a minimal barrier for protection against small amounts of fluid penetration. Fabric made in the USA from 65/35 Poly/Cotton.

## COMPLIANCE

- Level 1 fabrics/garments are:
  - Used for minimal risk situations
  - Provides a slight barrier to small amounts of fluid penetration
- AAMI PB70 Barrier Protection Level 1 verified by Vartest Laboratories
- Barrier fabrics can be laundered, sanitized and reused through a minimum of 25 industrial or home launderings
- Industrial wash @150 F per 8.1.3. of NFPA 2112-2018, no sterilization post washing
- This is not effective for the prevention or treatment of COVID-19
- The product has not been FDA cleared or approved
- The product has been authorized by FDA under an Emergency Use Authorization (EUA) for use as PPE in healthcare settings by HCP as they may help protect HCP and/or patients from the transfer of the SARS-CoV-2 virus in low or minimal risk level situations to prevent the spread of COVID-19
- This product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of medical devices, including alternative products used as medical devices, due to shortages during the COVID-19 outbreak, under Section 564(b)(1) of the Act, 21 USC 360bbb-3(b)(1) unless the authorization is terminated or revoked sooner

CODE	PRODUCT DESCRIPTION	PACKAGING	CASE COUNT	CASE DIM.(CM) LWH	CASE WT
406222	REUSABLE ISOLATION GOWN	3 GOWNS/BAG	16 BAGS	24 X 16 X 12	47 LBS

† Due to high demand assembled in Mexico or Nicaragua.

\*It is not intended to prevent or protect from any form of virus, bacteria, illness or disease (or otherwise). This product is being distributed according to the FDA's enforcement discretion policy applicable to non-surgical gowns and minimal-to-low barrier protection surgical apparel. The FDA does not object to the use of this gown in a healthcare setting provided the following:

- This gown may be used when FDA cleared gowns or apparel are unavailable.
- Not recommended for use in a surgical setting or where significant exposure to body fluids or other hazardous fluids may be expected.
- Not recommended for use in a clinical setting where Level 3 or Level 4 protection is warranted.
- Not recommended for use in the presence of high intensity heat source or flammable gas.